

Peptidyl-prolyl cis-trans isomerase A

Catalog No: #AP87316



Package Size: #AP87316-1 50ug #AP87316-2 100ug #AP87316-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Peptidyl-prolyl cis-trans isomerase A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PPlase A,Cyclophilin A,Cyclosporin A-binding protein,Rotamase A,p1B15,p31
Accession No.	Uniprot ID:P10111
Uniprot	P10111
GeneID	100360977;25518;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ppia

Note: This product is for in vitro research use only