

RING finger protein 39

Catalog No: #AP87319



Package Size: #AP87319-1 50ug #AP87319-2 100ug #AP87319-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RING finger protein 39
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	LTP-induced RING finger protein,Lirf
Accession No.	Uniprot ID:Q920M2
Uniprot	Q920M2
GeneID	171387;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Rnf39

Note: This product is for in vitro research use only