

POU domain, class 3, transcription factor 3

Catalog No: #AP87321



Package Size: #AP87321-1 50ug #AP87321-2 100ug #AP87321-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	POU domain, class 3, transcription factor 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Brain-specific homeobox/POU domain protein 1,Brain-1,Brn-1,Brn1,Rhs2
Accession No.	Uniprot ID:Q63262
Uniprot	Q63262
GeneID	192109;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pou3f3

Note: This product is for in vitro research use only