

Retinol-binding protein 4

Catalog No: #AP87322



Package Size: #AP87322-1 50ug #AP87322-2 100ug #AP87322-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Retinol-binding protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Plasma retinol-binding protein,PRBP
Accession No.	Uniprot ID:P41263
Uniprot	P41263
GeneID	396166;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RBP4

Note: This product is for in vitro research use only