

Resistin-like alpha

Catalog No: #AP87323



Package Size: #AP87323-1 50ug #AP87323-2 100ug #AP87323-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Resistin-like alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cysteine-rich secreted protein FIZZ1,RELMalpha,Fizz1
Accession No.	Uniprot ID:Q99P85
Uniprot	Q99P85
GeneID	81712;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Retnla

Note: This product is for in vitro research use only