

Prepronociceptin

Catalog No: #AP87324



Package Size: #AP87324-1 50ug #AP87324-2 100ug #AP87324-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Prepronociceptin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	N23K/N27K,Npnc1
Accession No.	Uniprot ID:Q64387
Uniprot	Q64387
GeneID	18155;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Pnoc

Note: This product is for in vitro research use only