

Retinol-binding protein 1

Catalog No: #AP87330



Package Size: #AP87330-1 50ug #AP87330-2 100ug #AP87330-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Retinol-binding protein 1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Bovine |
| Other Names | Cellular retinol-binding protein,CRBP,Cellular retinol-binding protein I,CRBP-I |
| Accession No. | Uniprot ID:P02694 |
| Uniprot | P02694 |
| GeneID | 537379; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RBP1

Note: This product is for in vitro research use only