

Ras-related protein Rab-8A

Catalog No: #AP87332



Package Size: #AP87332-1 50ug #AP87332-2 100ug #AP87332-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ras-related protein Rab-8A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Oncogene c-mel,MEL,RAB8
Accession No.	Uniprot ID:P61007
Uniprot	P61007
GeneID	403773;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RAB8A

Note: This product is for in vitro research use only