

RAF proto-oncogene serine/threonine-protein kinase

Catalog No: #AP87333



Package Size: #AP87333-1 50ug #AP87333-2 100ug #AP87333-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RAF proto-oncogene serine/threonine-protein kinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	C-MIL,C-RAF,MIL proto-oncogene serine/threonine-protein kinase,RAF-1,MIL
Accession No.	Uniprot ID:P05625
Uniprot	P05625
GeneID	396245;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RAF1

Note: This product is for in vitro research use only