

## Rho GTPase-activating protein 21

Catalog No: #AP87335



Package Size: #AP87335-1 50ug #AP87335-2 100ug #AP87335-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Rho GTPase-activating protein 21   |
| Brief Description  | Recombinant Protein  |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE   |
| Species Reactivity | Ms   |
| Other Names        | Rho GTPase-activating protein 10,Rho-type GTPase-activating protein 21,Arhgap10,Kiaa1424 |
| Accession No.      | Uniprot ID:Q6DFV3  |
| Uniprot            | Q6DFV3   |
| GeneID             | 71435;   |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)                                   |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Arhgap21

Note: This product is for in vitro research use only