

Pituitary tumor-transforming gene 1 protein-interacting protein

Catalog No: #AP87341

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP87341-1 50ug #AP87341-2 100ug #AP87341-3 1mg

Description

Product Name	Pituitary tumor-transforming gene 1 protein-interacting protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Pituitary tumor-transforming gene protein-binding factor,PBF
Accession No.	Uniprot ID:Q8R143
Uniprot	Q8R143
GeneID	108705;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pttg1ip

Note: This product is for in vitro research use only