

Sialin

Catalog No: #AP87346



Package Size: #AP87346-1 50ug #AP87346-2 100ug #AP87346-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sialin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	H(+)/nitrate cotransporter,H(+)/sialic acid cotransporter,AST,Solute carrier family 17 member 5,Vesicular H(+)/Aspartate-glutamate cotransporter
Accession No.	Uniprot ID:Q8BN82
Uniprot	Q8BN82
GeneID	235504;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Slc17a5

Note: This product is for in vitro research use only