

Proto-oncogene tyrosine-protein kinase receptor Ret

Catalog No: #AP87349



Package Size: #AP87349-1 50ug #AP87349-2 100ug #AP87349-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Proto-oncogene tyrosine-protein kinase receptor Ret
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cadherin family member 12,Proto-oncogene c-Ret,CDHF12,CDHR16,PTC,RET51
Accession No.	Uniprot ID:P07949
Uniprot	P07949
GeneID	5979;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RET

Note: This product is for in vitro research use only