

Adenosylhomocysteinase

Catalog No: #AP87352



Package Size: #AP87352-1 50ug #AP87352-2 100ug #AP87352-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Adenosylhomocysteinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	AdoHcyase,CUBP,Liver copper-binding protein,S-adenosyl-L-homocysteine hydrolase
Accession No.	Uniprot ID:P50247
Uniprot	P50247
GeneID	11615;269378;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ahcy

Note: This product is for in vitro research use only