

Secreted frizzled-related sequence protein 4

Catalog No: #AP87353



Package Size: #AP87353-1 50ug #AP87353-2 100ug #AP87353-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Secreted frizzled-related sequence protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	FRP-4
Accession No.	Uniprot ID:Q9Z1N6
Uniprot	Q9Z1N6
GeneID	20379;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Sfrp4

Note: This product is for in vitro research use only