

Serine/threonine-protein kinase Sgk2

Catalog No: #AP87354



Package Size: #AP87354-1 50ug #AP87354-2 100ug #AP87354-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/threonine-protein kinase Sgk2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Serum/glucocorticoid-regulated kinase 2
Accession No.	Uniprot ID:Q9HBY8
Uniprot	Q9HBY8
GeneID	10110;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SGK2

Note: This product is for in vitro research use only