

V(D)J recombination-activating protein 1

Catalog No: #AP87358



Package Size: #AP87358-1 50ug #AP87358-2 100ug #AP87358-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	V(D)J recombination-activating protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	RAG-1
Accession No.	Uniprot ID:Q867B5
Uniprot	Q867B5
GeneID	397506;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RAG1

Note: This product is for in vitro research use only