

RING-box protein 2

Catalog No: #AP87359



Package Size: #AP87359-1 50ug #AP87359-2 100ug #AP87359-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RING-box protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Rbx2,RING finger protein 7,Sensitive to apoptosis gene protein,Rbx2,Sag
Accession No.	Uniprot ID:Q9WTZ1
Uniprot	Q9WTZ1
GeneID	19823;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Rnf7

Note: This product is for in vitro research use only