

Ribose-5-phosphate isomerase

Catalog No: #AP87361



Package Size: #AP87361-1 50ug #AP87361-2 100ug #AP87361-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ribose-5-phosphate isomerase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Phosphoriboisomerase,Rpi
Accession No.	Uniprot ID:P47968
Uniprot	P47968
GeneID	19895;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Rpia

Note: This product is for in vitro research use only