

# Pituitary tumor-transforming gene 1 protein-interacting protein

Catalog No: #AP87363

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP87363-1 50ug #AP87363-2 100ug #AP87363-3 1mg

## Description

Product Name	Pituitary tumor-transforming gene 1 protein-interacting protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Pituitary tumor-transforming gene protein-binding factor,PBF,C21orf1,C21orf3
Accession No.	Uniprot ID:P53801
Uniprot	P53801
GeneID	754;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PTTG1IP

Note: This product is for in vitro research use only