

Sodium channel protein type 9 subunit alpha

Catalog No: #AP87365



Package Size: #AP87365-1 50ug #AP87365-2 100ug #AP87365-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sodium channel protein type 9 subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Neuroendocrine sodium channel,hNE-Na,Peripheral sodium channel 1,PN1,Sodium channel protein type IX subunit alpha,Voltage-gated sodium channel subunit alpha Nav1.7,NENA
Accession No.	Uniprot ID:Q15858
Uniprot	Q15858
GeneID	6335;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SCN9A

Note: This product is for in vitro research use only