

Reticulon-4

Catalog No: #AP87383



Package Size: #AP87383-1 50ug #AP87383-2 100ug #AP87383-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Reticulon-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Foocen, Neurite outgrowth inhibitor, Nogo protein, Neuroendocrine-specific protein, NSP, Neuroendocrine-specific protein C homolog, RTN-x, Reticulon-5, KIAA0886, NOGO, My043, SP1507
Accession No.	Uniprot ID:Q9NQC3
Uniprot	Q9NQC3
GeneID	57142;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: RTN4

Note: This product is for in vitro research use only