

Ras-related protein M-Ras

Catalog No: #AP87384



Package Size: #AP87384-1 50ug #AP87384-2 100ug #AP87384-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ras-related protein M-Ras
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Muscle and microspikes Ras,Ras-related protein R-Ras3,X-Ras,Xras
Accession No.	Uniprot ID:O08989
Uniprot	O08989
GeneID	17532;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mras

Note: This product is for in vitro research use only