

Peroxiredoxin-2

Catalog No: #AP87386



Package Size: #AP87386-1 50ug #AP87386-2 100ug #AP87386-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Peroxiredoxin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Thiol-specific antioxidant protein, TSA, Thioredoxin peroxidase 1, Thioredoxin-dependent peroxide reductase 1, Tdp1, Tpx
Accession No.	Uniprot ID:Q61171
Uniprot	Q61171
GeneID	21672;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Prdx2

Note: This product is for in vitro research use only