

Peroxisome proliferator-activated receptor gamma coactivator 1-alpha

Catalog No: #AP87387

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP87387-1 50ug #AP87387-2 100ug #AP87387-3 1mg

Description

Product Name	Peroxisome proliferator-activated receptor gamma coactivator 1-alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PGC-1-alpha,Ligand effect modulator 6,LEM6,PGC1,PGC1A,PPARGC1
Accession No.	Uniprot ID:Q9UBK2
Uniprot	Q9UBK2
GeneID	10891;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PPARGC1A

Note: This product is for in vitro research use only