

## Rho-associated protein kinase 2

Catalog No: #AP87396



Package Size: #AP87396-1 50ug #AP87396-2 100ug #AP87396-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rho-associated protein kinase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Rho-associated, coiled-coil-containing protein kinase 2,Rho-associated, coiled-coil-containing protein kinase II,ROCK-II,p164 ROCK-2
Accession No.	Uniprot ID:Q28021
Uniprot	Q28021
GeneID	282041;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ROCK2

Note: This product is for in vitro research use only