

RuvB-like 1

Catalog No: #AP87402

Package Size: #AP87402-1 50ug #AP87402-2 100ug #AP87402-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RuvB-like 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	49 kDa TATA box-binding protein-interacting protein,49 kDa TBP-interacting protein,54 kDa erythrocyte cytosolic protein,ECP-54,INO80 complex subunit H,Nuclear matrix protein 238,NMP 238,Pontin 52,TIP49a,TIP60-associated protein 54-alpha,TAP54-alpha,INO80H,NMP238,TIP49,TIP49A
Accession No.	Uniprot ID:Q9Y265
Uniprot	Q9Y265
GeneID	8607;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RUVBL1

Note: This product is for in vitro research use only