

Parathyroid hormone/parathyroid hormone-related peptide receptor

Catalog No: #AP87405

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP87405-1 50ug #AP87405-2 100ug #AP87405-3 1mg

Description

Product Name	Parathyroid hormone/parathyroid hormone-related peptide receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	PTH/PTHrP type I receptor,PTH/PTHr receptor,Parathyroid hormone 1 receptor,PTH1 receptor,PTHR,PTHR1
Accession No.	Uniprot ID:P50133
Uniprot	P50133
GeneID	397675;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PTH1R

Note: This product is for in vitro research use only