

Retinoblastoma-like protein 2

Catalog No: #AP87406



Package Size: #AP87406-1 50ug #AP87406-2 100ug #AP87406-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Retinoblastoma-like protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	130 kDa retinoblastoma-associated protein,p130,PPAR-alpha-interacting complex protein 128,PRIC128,Retinoblastoma-related protein 2,RBR-2,pRb2
Accession No.	Uniprot ID:O55081
Uniprot	O55081
GeneID	81758;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Rb12

Note: This product is for in vitro research use only