

Meiotic recombination protein REC8 homolog

Catalog No: #AP87407



Package Size: #AP87407-1 50ug #AP87407-2 100ug #AP87407-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Meiotic recombination protein REC8 homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cohesin Rec8p,Rec8L1
Accession No.	Uniprot ID:Q6AYJ4
Uniprot	Q6AYJ4
GeneID	290227;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Rec8

Note: This product is for in vitro research use only