

Double-stranded RNA-specific editase 1

Catalog No: #AP87408



Package Size: #AP87408-1 50ug #AP87408-2 100ug #AP87408-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Double-stranded RNA-specific editase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	RNA-editing deaminase 1, RNA-editing enzyme 1, dsRNA adenosine deaminase, Adar2, Red1
Accession No.	Uniprot ID:Q91ZS8
Uniprot	Q91ZS8
GeneID	110532;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Adarb1

Note: This product is for in vitro research use only