

Poliovirus receptor-related protein 1

Catalog No: #AP87410



Package Size: #AP87410-1 50ug #AP87410-2 100ug #AP87410-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Poliovirus receptor-related protein 1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | Herpes virus entry mediator C, Herpesvirus entry mediator C, HveC, Nectin-1, Hvec, Prr1, CD111 |
| Accession No. | Uniprot ID:Q9JKF6 |
| Uniprot | Q9JKF6 |
| GeneID | 58235; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Pvr1

Note: This product is for in vitro research use only