

Secretory carrier-associated membrane protein 1

Catalog No: #AP87414



Package Size: #AP87414-1 50ug #AP87414-2 100ug #AP87414-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Secretory carrier-associated membrane protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Secretory carrier membrane protein 1,SCAMP 37,Scamp
Accession No.	Uniprot ID:P56603
Uniprot	P56603
GeneID	29521;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Scamp1

Note: This product is for in vitro research use only