

# Reactive oxygen species modulator 1

Catalog No: #AP87415



Package Size: #AP87415-1 50ug #AP87415-2 100ug #AP87415-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                        |
|--------------------|--------------------------------------------------------|
| Product Name       | Reactive oxygen species modulator 1                    |
| Brief Description  | Recombinant Protein                                    |
| Host Species       | E.coli                                                 |
| Purification       | Greater than 90% by SDS-PAGE                           |
| Species Reactivity | Bovine                                                 |
| Other Names        | ROS modulator 1,Protein MGR2 homolog                   |
| Accession No.      | Uniprot ID:Q3SZV8                                      |
| Uniprot            | Q3SZV8                                                 |
| GeneID             | 618530;                                                |
| Tag Info           | His                                                    |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.) |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ROMO1

Note: This product is for in vitro research use only