

# Polypyrimidine tract-binding protein 1

Catalog No: #AP87416



Package Size: #AP87416-1 50ug #AP87416-2 100ug #AP87416-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Polypyrimidine tract-binding protein 1                    |
| Brief Description  | Recombinant Protein                                       |
| Host Species       | E.coli  |
| Purification       | Greater than 90% by SDS-PAGE                              |
| Species Reactivity | Ms  |
| Other Names        | PTB,Heterogeneous nuclear ribonucleoprotein I,hnRNP I,Ptb |
| Accession No.      | Uniprot ID:P17225   |
| Uniprot            | P17225  |
| Tag Info           | His   |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)    |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ptbp1

Note: This product is for in vitro research use only