

Tyrosine-protein phosphatase non-receptor type 1

Catalog No: #AP87417



Package Size: #AP87417-1 50ug #AP87417-2 100ug #AP87417-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tyrosine-protein phosphatase non-receptor type 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	CPTP1,Protein-tyrosine phosphatase 1B,PTP-1B,PTP1B
Accession No.	Uniprot ID:O13016
Uniprot	O13016
GeneID	395688;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PTPN1

Note: This product is for in vitro research use only