

# Major prion protein

Catalog No: #AP87421



Package Size: #AP87421-1 50ug #AP87421-2 100ug #AP87421-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Major prion protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PrP,CD230,Pm,Prp
Accession No.	Uniprot ID:P13852
Uniprot	P13852
GeneID	24686;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Prnp

Note: This product is for in vitro research use only