

Pregnancy zone protein

Catalog No: #AP87423



Package Size: #AP87423-1 50ug #AP87423-2 100ug #AP87423-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pregnancy zone protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	C3 and PZP-like alpha-2-macroglobulin domain-containing protein 6,CPAMD6
Accession No.	Uniprot ID:P20742
Uniprot	P20742
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PZP

Note: This product is for in vitro research use only