

Resistin-like alpha

Catalog No: #AP87427



Package Size: #AP87427-1 50ug #AP87427-2 100ug #AP87427-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Resistin-like alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cysteine-rich secreted protein A12-gamma, Cysteine-rich secreted protein FIZZ1, Hypoxia-induced mitogenic factor, Parasite-induced macrophage novel gene 1 protein, RELMalpha, Fizz1, Himf, Pmng1
Accession No.	Uniprot ID:Q9EP95
Uniprot	Q9EP95
GeneID	57262;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Retnla

Note: This product is for in vitro research use only