

RNA exonuclease 1 homolog

Catalog No: #AP87428



Package Size: #AP87428-1 50ug #AP87428-2 100ug #AP87428-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | RNA exonuclease 1 homolog |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | Elongin-A-binding protein 1,EloA-BP1,Transcription elongation factor B polypeptide 3-binding protein 1,ELOABP1,KIAA1138,TCEB3BP1 |
| Accession No. | Uniprot ID:Q8N1G1 |
| Uniprot | Q8N1G1 |
| GeneID | 57455; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:REXO1

Note: This product is for in vitro research use only