

## Splicing factor, proline- and glutamine-rich

Catalog No: #AP87432



Package Size: #AP87432-1 50ug #AP87432-2 100ug #AP87432-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Splicing factor, proline- and glutamine-rich
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	100 kDa DNA-pairing protein,hPOMp100,DNA-binding p52/p100 complex, 100 kDa subunit,Polypyrimidine tract-binding protein-associated-splicing factor,PSF,PSF
Accession No.	Uniprot ID:P23246
Uniprot	P23246
GeneID	6421;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SFPQ

Note: This product is for in vitro research use only