

Hemojuvelin

Catalog No: #AP87434



Package Size: #AP87434-1 50ug #AP87434-2 100ug #AP87434-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hemojuvelin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Hemochromatosis type 2 protein homolog,RGM domain family member C,Rgmc
Accession No.	Uniprot ID:Q7TQ32
Uniprot	Q7TQ32
GeneID	69585;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hfe2

Note: This product is for in vitro research use only