

Renin-1

Catalog No: #AP87435



Package Size: #AP87435-1 50ug #AP87435-2 100ug #AP87435-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Renin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Angiotensinogenase,Kidney renin,Ren,Ren-1
Accession No.	Uniprot ID:P06281
Uniprot	P06281
GeneID	19701;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ren1

Note: This product is for in vitro research use only