

Sodium channel protein type 8 subunit alpha

Catalog No: #AP87444



Package Size: #AP87444-1 50ug #AP87444-2 100ug #AP87444-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sodium channel protein type 8 subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Sodium channel protein type VIII subunit alpha, Voltage-gated sodium channel subunit alpha Nav1.6, MED
Accession No.	Uniprot ID:Q9UQD0
Uniprot	Q9UQD0
GeneID	6334;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SCN8A

Note: This product is for in vitro research use only