

## Proteasome subunit beta type-9

Catalog No: #AP87446



Package Size: #AP87446-1 50ug #AP87446-2 100ug #AP87446-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Proteasome subunit beta type-9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	LMP-2d,Low molecular mass protein 2,Macropain chain 7,Multicatalytic endopeptidase complex chain 7,Proteasome chain 7,Proteasome subunit beta-1i,Really interesting new gene 12 protein,Lmp2,Ring12
Accession No.	Uniprot ID:P28076
Uniprot	P28076
GeneID	16912;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Psmb9

Note: This product is for in vitro research use only