

Regulator of G-protein signaling 9

Catalog No: #AP87447



Package Size: #AP87447-1 50ug #AP87447-2 100ug #AP87447-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Regulator of G-protein signaling 9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	RGS9
Accession No.	Uniprot ID:O54828
Uniprot	O54828
GeneID	19739;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Rgs9

Note: This product is for in vitro research use only