

Rho-associated protein kinase 2

Catalog No: #AP87450



Package Size: #AP87450-1 50ug #AP87450-2 100ug #AP87450-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rho-associated protein kinase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Rho-associated, coiled-coil-containing protein kinase 2,Rho-associated, coiled-coil-containing protein kinase II,ROCK-II,RhoA-binding kinase 2,p150 ROK-alpha,ROKalpha,p164 ROCK-2
Accession No.	Uniprot ID:Q62868
Uniprot	Q62868
GeneID	25537;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Rock2

Note: This product is for in vitro research use only