

Protein S100-A4

Catalog No: #AP87453



Package Size: #AP87453-1 50ug #AP87453-2 100ug #AP87453-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein S100-A4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Metastasin, Nerve growth factor-induced protein 42A, P9K, Placental calcium-binding protein, S100 calcium-binding protein A4
Accession No.	Uniprot ID:P05942
Uniprot	P05942
GeneID	24615;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: S100a4

Note: This product is for in vitro research use only