

Receptor-interacting serine/threonine-protein kinase 2

Catalog No: #AP87457



Package Size: #AP87457-1 50ug #AP87457-2 100ug #AP87457-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Receptor-interacting serine/threonine-protein kinase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CARD-containing interleukin-1 beta-converting enzyme-associated kinase,CARD-containing IL-1 beta ICE-kinase,RIP-like-interacting CLARP kinase,Receptor-interacting protein 2,RIP-2,Tyrosine-protein kinase RIPK2,CARDIAK,RICK,RIP2,UNQ277/PRO314/PRO34092
Accession No.	Uniprot ID:O43353
Uniprot	O43353
GeneID	8767;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RIPK2

Note: This product is for in vitro research use only