

Retinoic acid receptor RXR-beta

Catalog No: #AP87458



Package Size: #AP87458-1 50ug #AP87458-2 100ug #AP87458-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Retinoic acid receptor RXR-beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Nuclear receptor subfamily 2 group B member 2, Retinoid X receptor beta, NR2B2
Accession No.	Uniprot ID:Q5TJF7
Uniprot	Q5TJF7
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: RXRB

Note: This product is for in vitro research use only